

MOLTOX[®]

Molecular Toxicology, Inc.

Phenotype Test Kit Instruction Manual

31-600



INTRODUCTION

Phenotype tests are performed to verify the genotype of *S. typhimurium* and *E. coli* strains utilized in the reverse bacterial mutation test. This is a requirement of the OECD 471 guideline. As opposed to traditional tests that require multiple plates per strain and specialized equipment, this kit allows for use of one plate per strain, can be used with both *S. typhimurium* and *E. coli* tester strains, and utilizes Mitomycin C discs in place of UV light exposure. Results are available the next day.

Principle

Wells 1 – 4 contain Minimal Glucose Agar (MGA) with excess histidine, biotin, and tryptophan. The strain culture is added to each well and a specific phenotype test disc applied to each. Discs and the characteristic they detect are as follows:

Designation	Disc type	Detects	Genotype
C	Crystal Violet	Sensitivity to crystal violet	<i>rfa</i> deep mutation
A	Ampicillin	Ampicillin resistance	+/- pKM101 R-factor plasmid
T	Tetracycline	Tetracycline resistance	+/- pAQ1 plasmid
M	Mitomycin C	Sensitivity to mitomycin C ¹	<i>uvrA</i> and <i>uvrB</i> repair deficiency ¹

¹ Traditionally, *uvrA* and *uvrB* repair deficiency is tested via exposure to UV light, a cross-linking agent. Use of mitomycin C, also a cross-linking agent, is used in place of UV exposure.

Wells 5 and 6 contain MGA with excess histidine and biotin and MGA with excess tryptophan and biotin, respectively. These wells test for tryptophan and histidine growth requirements.

Additional items you will need to complete the assay

Top agar, non-supplemented (≥ 8 mls)	Water bath, 45°C - 47°C
Sharpie or wax pencil	Test tubes, 13 x 100 mm, sterile (2)
Micropipettes, 100 µL - 500 µL	Test tube rack
Micropipette tips, 100 µL - 1000 µL, sterile	Incubator, 37°C
Pipet-aid or rubber bulb	Vortex mixer
Pipets, 5 or 10 mL, sterile	Forceps, small, sterile
Microwave or boiling water bath	

Procedure

Proceed as below for each strain using sterile technique

Use aseptic technique throughout the procedure

1. Using a sharpie or wax pencil, label plate(s) on the side or lid with the strain designation(s).
2. Loosen the top agar container cap then melt the top agar in a boiling water bath or microwave oven. Examine the melted agar carefully. If any opalescence persists, continue heating until a perfectly clear solution is obtained. After melting, place the top agar bottle into a 45°C - 47°C water bath. Allow at least 45 minutes for temperature equilibration.
3. Add 4 mls molten top agar to a 13 x 100 mm test tube. Add 100 µL of bacterial suspension culture to the top agar and mix.
4. Dispense 0.5 mL of the mix into each well of the 6-well phenotype plate ensuring it covers the surface of each well. Leave the plate on a level surface for the agar to gel.
5. Using forceps, place one C, A, T and M (crystal violet, ampicillin, tetracycline and mitomycin C) disc near the center of wells 1 to 4 respectively.
6. Invert the plate(s) and incubate at 37°C for 16 - 24 hours.
7. Hold the plate 3 or 4 cm above a printed page then estimate the diameter of the clear zone of inhibition in wells 1 to 4 using a ruler. Record density of growth in wells 5 and 6.

Expected results

Strain	Zone of inhibition (mm)				Growth	
	C	A	T	M	<i>trp</i> ⁻	<i>his</i> ⁻
TA1535, TA1537, TA1538	x	x	x	x	++	-
TA97a, TA98, TA100	x		x	x	++	-
TA102	x				++	-
WP2		x			++	++
WP2 <i>uvrA</i>		x	x	x	-/+	++
WP2 pKM101			x		-/+	++
WP2 <i>uvrA</i> pKM101			x	x	-/+	++
	x	zone > 10 mm	-	no growth	-/+	no or weak growth
			+	weak growth	++	strong growth

Version

28 April 2026

MOLTOX[®]
Molecular Toxicology, Inc.

157 Industrial Park Drive Boone, NC 28607 828.264.9099 Toll Free: 800.536.7232 Fax: 828.264.0103

www.MOLTOX.com sales@moltox.com