


Genotype Confirmation Testing of *E. coli* WP2 Ames Strains

		Assay Setup	Assay Performance	Expected Results
Genotype	Test Media	Day 0	Day 1	Day 2
<i>trpE</i>	<p>MGA¹</p> <p>&</p> <p>MGA¹ w/ 10 - 50 µg/mL tryptophan</p>	Inoculate cultures ² (15 - 18 h prior to assay performance)	<p>Grow cultures to 1 - 2 x 10⁹ cfu/ml²</p> <p>Using a 10 µl inoculating loop³ or a swab^{3,4}, streak culture(s) across MGA(s)</p>  <p>Incubate <i>overnight</i>, 37°C</p>	No growth on MGA; growth on MGA w/tryptophan
<i>pKM101</i>⁵	MGA ¹ w/ 10 - 50 µg/mL tryptophan & 24 µg/mL ampicillin	Prepare phenotype plates (if applicable)	<p>Incubate <i>overnight</i>, 37°C</p>	Growth of <i>pKM101</i> - containing strains only ²
<i>uvrA</i>^{6,7}	MGA ¹ w/ 10 - 50 µg/mL tryptophan ⁸			No growth of <i>uvrA</i> -containing strains on +UV side ^{2, 6, 7}

¹ Minimal Glucose Agar. Reference Mortelmans and Zeiger (2000), The Ames *Salmonella*/microsome mutagenicity assay, *Mut Res.* 455, pp 29-70.

² Reference MOLTOX[®] Technical Bulletin *Lyophilized Culture Discs*.

³ Sterile

⁴ Press swab against the side of the growth vessel to remove excess broth prior to inoculating MGA plate. This helps prevent nutrient carryover to the MGA plate which can cause erroneous results.

⁵ Streak a *pKM101*-containing and a *pKM101*-free strain on the same plate for growth comparison after incubation. Reference MOLTOX[®] Technical Bulletin *Lyophilized Culture Discs*.

⁶ Streak a WP2 and a WP2 *uvrA* culture(s) across plate². With dish lid off, cover ½ the plate and expose the uncovered side to a 15W UV lamp at a distance of 33 cm. Expose non-*pKM101* strains for 6 sec, *pKM101*-containing strains for 8 sec.

⁷ Alternate *uvrA* test: After streaking, apply a filter paper disc containing 0.2 µg Mitomycin C to the center of the streak. *uvrA*-containing strains will exhibit a halo around the disc after incubation.

⁸ Nutrient Agar prepared with Oxoid Broth No.2 can also be used.